

PerOssal® 1x50 pcs.

1 unit with 50 pellets is sufficient for filling of 10 cc


Gentamicin
Rifampicin

The treating physician is responsible for the decision regarding the type and quantity of the corresponding antibiotic.


Combined Application

Preparation


Dissolve 2 ampoules Gentamicin 120 mg and 300 mg Rifampicin in 5 ml



Dissolve 2 ampoules Gentamicin 120 mg and 600 mg Rifampicin in 10 ml



Dissolve 4 ampoules Gentamicin 120 mg and 600 mg Rifampicin in 10 ml




Dissolve 4 ampoules Gentamicin 120 mg and 600 mg Rifampicin in 5 ml and 500 mg Rifampicin in 10 ml




Drenching


240 mg Gentamicin 240 mg Rifampicin 2 units PerOssal® (20 cc)




240 mg Gentamicin 480 mg Rifampicin 3 units PerOssal® (30 cc)



480 mg Gentamicin 480 mg Rifampicin 4 units PerOssal® (40 cc)



480 mg Gentamicin 720 mg Rifampicin 5 units PerOssal® (50 cc)



PerOssal® 1x6 & 2x6 pcs.

1 tube with 6 pellets is sufficient for filling of 1.2 cc

Gentamicin
or
Tobramycin

The dosage recommendation refers to the in vitro results, as well as to the results obtained within the scope of pilot testing. The treating physician is responsible for the final decision regarding the type and quantity of the corresponding antibiotic.

Single Application

Preparation

1 ampoule Gentamicin or Tobramycin 40 mg/ml




1 ampoule Gentamicin or Tobramycin 40 mg/ml




Drenching

20 mg Gentamicin or Tobramycin 1x6 pellets PerOssal® (1.2 cc)




40 mg Gentamicin or Tobramycin 2x6 pellets PerOssal® (2.4 cc)




Preparation

Dissolve 500 mg Vancomycin as a dry substance in 5 ml




Dissolve 500 mg Vancomycin as a dry substance in 5 ml




Drenching

50 mg Vancomycin 1x6 pellets PerOssal® (1.2 cc)



100 mg Vancomycin 2x6 pellets PerOssal® (2.4 cc)



aap Implantate AG
Kronenweg 10 • 42699 Solingen • Germany
Phone: +49 201 75013-0 • Fax: +49 201 75013-222
customerservice@aap.de
www.aap.de

responsible manufacturer:
aap Biomaterials GmbH
Kronenweg 10 • 42699 Solingen • Germany
Phone: +49 201 75013-0 • Fax: +49 201 75013-100
biomaterials@aap.de
www.aap.de



Dosage Recommendation

PerOssal®
Resorbable Bone Matrix

PerOssal® 1x50 pcs.

1 unit with 50 pellets is sufficient for filling of 10 cc

The dosage recommendation refers to the in vitro results, as well as to the results obtained within the scope of their testing. The treating physician is responsible for the final decision regarding the type and quantity of the corresponding antibiotic.

Single Application

Gentamicin or Tobramycin

Preparation

2 ampoules Gentamicin
or Tobramycin 80 mg

2x2 ml

Drenching

► 160 mg Gentamicin or Tobramycin ► 1 unit PerOssal® (10 cc)

4 ml

4 ampoules Gentamicin
or Tobramycin 80 mg

4x2 ml

► 320 mg Gentamicin or Tobramycin ► 2 units PerOssal® (20 cc)

4 ml 4 ml

6 ampoules Gentamicin
or Tobramycin 80 mg

6x2 ml

► 480 mg Gentamicin or Tobramycin ► 3 units PerOssal® (30 cc)

4 ml 4 ml 4 ml

6 ampoules Gentamicin
or Tobramycin 80 mg

6x2 ml

► 480 mg Gentamicin or Tobramycin ► 4 units PerOssal® (40 cc)

4 ml 4 ml 4 ml empty

6 ampoules Gentamicin
or Tobramycin 80 mg

6x2 ml

► 480 mg Gentamicin or Tobramycin ► 5 units PerOssal® (50 cc)

4 ml 4 ml 4 ml empty

PerOssal® 1x50 pcs.

1 unit with 50 pellets is sufficient for filling of 10 cc

The dosage recommendation refers to the in vitro results, as well as to the results obtained within the scope of their testing. The treating physician is responsible for the final decision regarding the type and quantity of the corresponding antibiotic.

Single Application

Vancomycin

Preparation

Dissolve 500 mg
Vancomycin as a
dry substance
in 5 ml

Drenching

► 400 mg Vancomycin ► 1 unit PerOssal® (10 cc)

4 ml

Dissolve 500 mg
Vancomycin as a
dry substance
in 10 ml

► 400 mg Vancomycin ► 2 units PerOssal® (20 cc)

4 ml 4 ml

Dissolve 500 mg
Vancomycin as a
dry substance
in 10 ml

► 500 mg Vancomycin ► 3 units PerOssal® (30 cc)

4 ml 4 ml 2 ml

Dissolve 2 units
with 500 mg
Vancomycin
each as a dry
substance
in 10 ml

► 800 mg Vancomycin ► 4 units PerOssal® (40 cc)

4 ml 4 ml 2x2 ml

Dissolve 2 units
with 500 mg
Vancomycin
each as a dry
substance
in 10 ml

► 1000 mg Vancomycin ► 5 units PerOssal® (50 cc)

4 ml 4 ml 2x2 ml

PerOssal® 1x50 pcs.

1 unit with 50 pellets is sufficient for filling of 10 cc

The dosage recommendation refers to the in vitro results, as well as to the results obtained within the scope of their testing. The treating physician is responsible for the final decision regarding the type and quantity of the corresponding antibiotic.

Combined Application

Preparation

Dissolve 500 mg Vancomycin
in 5 ml and
300 mg Rifampicin
in 5 ml

Drenching

► 400 mg Vancomycin ► 240 mg Rifampicin ► 2 units PerOssal® (20 cc)

4 ml

4 ml

Dissolve 500 mg Vancomycin
in 10 ml and
600 mg Rifampicin
in 10 ml

► 400 mg Vancomycin ► 480 mg Rifampicin ► 3 units PerOssal® (30 cc)

4 ml

4 ml 4 ml

Dissolve 1000 mg Vancomycin
in 10 ml and
600 mg Rifampicin
in 10 ml

► 800 mg Vancomycin ► 480 mg Rifampicin ► 4 units PerOssal® (40 cc)

4 ml

4 ml 4 ml

Dissolve 1000 mg Vancomycin
in 10 ml and
800 mg Rifampicin
in 10 ml

► 800 mg Vancomycin ► 720 mg Rifampicin ► 5 units PerOssal® (50 cc)

4 ml

4 ml 4 ml 2x2 ml